

Tabel de conversie ochelar – lentilă de contact

Minus

Ochelar	Lentilă contact	Ochelar	Lentilă contact
Vertex 12 mm	Vertex 0 mm	Vertex 12 mm	Vertex 0 mm
-3.50	-3.36	-20.00	-16.13
-3.75	-3.59	-20.50	-16.45
-4.00	-3.82	-21.00	-16.77
-4.25	-4.04	-21.50	-17.09
-4.50	-4.27	-22.00	-17.41
-4.75	-4.49	-22.50	-17.72
-5.00	-4.72	-23.00	-18.03
-5.25	-4.94	-23.50	-18.33
-5.50	-5.16	-24.00	-18.63
-5.75	-5.38	-24.50	-18.93
-6.00	-5.60	-25.00	-19.23
-6.25	-5.81	-25.50	-19.53
-6.50	-6.03	-26.00	-19.82
-6.75	-6.24	-26.50	-20.11
-7.00	-6.46	-27.00	-20.39
-7.25	-6.67	-27.50	-20.68
-7.50	-6.88	-28.00	-20.96
-7.75	-7.09	-28.50	-21.24
-8.00	-7.30	-29.00	-21.51
-8.25	-7.51	-29.50	-21.79
-8.50	-7.71	-30.00	-22.06
-8.75	-7.92	-30.50	-22.33
-9.00	-8.12	-31.00	-22.59
-9.25	-8.33	-31.50	-22.86
-9.50	-8.53	-32.00	-23.12
-9.75	-8.73	-32.50	-23.38
-10.00	-8.93	-33.00	-23.64
-10.50	-9.33	-33.50	-23.89
-11.00	-9.72	-34.00	-24.15
-11.50	-10.11	-34.50	-24.40
-12.00	-10.49	-35.00	-24.65
-12.50	-10.87	-35.50	-24.89
-13.00	-11.25	-36.00	-25.14
-13.50	-11.62		
-14.00	-11.99		
-14.50	-12.35		
-15.00	-12.71		
-15.50	-13.07		
-16.00	-13.42		
-16.50	-13.77		
-17.00	-14.12		
-17.50	-14.46		
-18.00	-14.80		
-18.50	-15.14		
-19.00	-15.47		
-19.50	-15.80		

Plus

Ochelar	Lentilă contact
Vertex 12 mm	Vertex 0 mm
3.50	3.65
3.75	3.93
4.00	4.20
4.25	4.48
4.50	4.76
4.75	5.04
5.00	5.32
5.25	5.60
5.50	5.89
5.75	6.18
6.00	6.47
6.25	6.76
6.50	7.05
6.75	7.34
7.00	7.64
7.25	7.94
7.50	8.24
7.75	8.54
8.00	8.85
8.25	9.16
8.50	9.47
8.75	9.78
9.00	10.09
9.25	10.40
9.50	10.72
9.75	11.04
10.00	11.36
10.50	12.01
11.00	12.67
11.50	13.34
12.00	14.02
12.50	14.71
13.00	15.40
13.50	16.11
14.00	16.83
14.50	17.55
15.00	18.29
15.50	19.04
16.00	19.80
16.50	20.57
17.00	21.36
17.50	22.15
18.00	22.96
18.50	23.78
19.00	24.61
19.50	25.46